


The Mysteries of the Injured Brain: Can they be solved with network physiology?


Dick Moberg
Moberg ICU Solutions

Traumatic Brain Injury (U.S.)



- Incidence:
 - 2.5 million/year
 - 250,000 hospitalized
 - 50,000 deaths
- Cost of TBI:
 - 76 billion dollars
 - 5 million people living with a TBI-related disability
- Conclusion
 - Small improvement in management/outcome can have a high impact on quality of life and cost of disease


Connecting the Dots



Steve Jobs Commencement Speech Stanford - 2005




Two Parallel Stories




1970s →

Traumatic Brain Injury

Brain Monitoring



CNS Injury – 1970s



Injury



Brain

Damage is done, nothing we can do


Spinal Cord

Impact + Secondary Injury

Spine Injury – 1970s




Jewell Osterholm



The pathophysiological response to spinal cord injury

Jewell L. Osterholm, *Journal of Neurosurgery*, January 1974 / Vol. 40 / No. 1 : Pages 3-33

Brain Injury – 1970s



Monitoring

Very crude metrics

- GCS (1974)
- ICP (Richmond bolt)
- Pupils


Imaging

- None

Secondary Injury

- “None”
- There is nothing we can do for these patients

Two Parallel Stories




1970s →

Secondary Injury in Spinal Cord

Traumatic Brain Injury

Brain Monitoring



Brain Monitoring 1970s – 1980s




1970

1980s

1990

2000

Nicolet
Med-80



Jefferson
Hospital

Neurotrac




Interspec


Brain Monitoring 1970s – 1980s



EEG
Cortical activity



Brain Monitoring 1970s – 1980s



Spectral Edge Frequency

1st EEG Metric

Used for Monitoring


Thiopental Induction

40 y.o. W. Male


Fp1 - T4 Lead

80 μ V P-P

16 sec/line



Brain Monitoring 1990s




1989



Experimental Device

NIH/NINDS Funding



1994



Neurotrac II


Continuous Monitoring

EEG + EP + TCD

EEG Metrics

- Power Bands (Absolute, Percent)
- Mean Dominant Frequency
- Median Power Frequency
- Peak Power Frequency
- Spectral Edge Frequency

Two Parallel Stories



1970s



Secondary Injury in Spinal Cord


Traumatic Brain Injury

Brain Monitoring


Single Metrics

- ICP
- EEG
- GCS

Remaining Work



CNS Injury – 1980s – 1990s



Injury




Brain

Impact + Secondary Injury

Spinal Cord

Impact + Secondary Injury

Brain Injury – 1990s



Monitoring

Very crude metrics

- GCS (1974)
- ICP (Richmond bolt)
- Pupils


Imaging

- CT (early 1980s)

Secondary Injury

- Apoptosis
- Many trials started
 - **All failed**

Why All the Failures?



- Complexity of the brain
- Too much variability in TBI patients

- Narayan et. al. - 2002

Two Parallel Stories



1970s



Secondary Injury in Spinal Cord

Secondary Injury In Brain but No Successful Trials


Traumatic Brain Injury

Brain Monitoring


Single Metrics

- ICP
- EEG
- GCS

Remaining Work



Multimodal Monitoring



Gert Pfurtscheller



Graz




Karin





Gerhard Litscher




Brain Roadmap – Dr. Gerhard Litscher



Brain Roadmap – System Connections




Multimodal Monitoring - 1996



Multimodal Monitoring

- Grant from NIH
- Developed computer architecture for gathering data from multiple products

Brain Monitoring 1970s – 2000s




1970

1980s

1990

2000

Nicolet
Med-80



Jefferson
Hospital

Neurotrac



Interspec

Neurotrac II




Moberg
Medical

CNS Monitor



Moberg
Research

The Product



The Component Neuromonitoring System (CNS)




A System of Components for Managing Data in
Neurocritical Care




The Smart Neuro ICU

Traditional Vital Signs Monitoring



Vital Signs Monitoring

- ECG
- HR
- ABP
- Temp
- CO
- RR
- SpO2




Challenge: Current Neuromonitoring



Vital Signs Monitoring


ECG
HR
ABP
Temp
CO
RR
SpO2



Neuromonitoring: Today




Solution




The vital signs monitor for the brain.



30+ devices
140+ measurements



Collect Data from Brain/Body Systems



Metabolism



Brain Oxygen



EEG



Perfusion



TCD



Temperature



Vital Signs



Pupils



Sedation



ICP



Video



Hemodynamics




CO2



Device - Pupilometer




Data Recording and Display




30+ Devices



Time-Synchronized Data



Typical CNS Installation




■ **MOBERG**
CNS Monitor

Vital Signs
Monitor



Dr. Michael Stiefel, Chair of Neurosurgery
Westchester Medical Center, Valhalla, NY

Brain System Interactions





Courtesy: Howard Yonas, MD, Univ of New Mexico

Missing Data – Wrong Treatment




Missing Data – Wrong Treatment



Courtesy of Columbia University

Two Parallel Stories



1970s



Secondary Injury in Spinal Cord

Secondary Injury In Brain but No Successful Trials


Traumatic Brain Injury

Brain Monitoring


Single Metrics
- ICP
- EEG
- GCS

Multimodal Monitoring

Remaining Work




TRACK-TBI and CENTER TBI




TRACK-TBI (large multi-center TBI trial) will produce advances in the **classification** of TBI patients as well as **outcome assessment**. It ignores the variability in patient management from site to site and from nursing shift to nursing shift.

Classification



Biomarkers



Advanced
Imaging

Genomics

Management Assessment



The CNS Monitor guides
the management of brain
injured patients.


Outcome Assessment

Psychometric
Testing


TED

Traumatic
Endpoints
Development


2014 TRACK-TBI and CENTER-TBI




Two large multi-center trials focused on characterizing TBI patients such that therapeutic measures can be targeted to the appropriate subgroup.




11 sites
3,000 subjects



80 sites, 21 countries
5,400 subjects




Big Data for a Complex Problem





- Lots of data being recorded
- Promise of progress with TBI
- Still a very complex problem to solve
- But how can we accelerate progress
 - How can we get innovative ideas to the clinics quickly


Two Parallel Stories



1970s



Electronic Medical Records




Larry Weed
(1923-2017)

Father of the Electronic Medical Record

Medical records that guide and teach
NEJM. **278** (11): 593–600.

Met in 1980. An “isolated dot”

Electronic Medical Records





Visit Date	POV Narrative	ICD	ICD N	ICD N
07/22/2005	Fractured lens	3999	Uncon	
07/22/2005	TYPE 2 DIABETES MELLITUS	250.00	Uncon	
07/19/2005	HYPERTENSION	401.9	Uncon	
07/19/2005	Fractured lens	3999	Uncon	
06/29/2005	TYPE 2 DIABETES MELLITUS	250.00	Uncon	
06/29/2005	HYPERTENSION	401.9	Uncon	
05/16/2005	TYPE 2 DIABETES MELLITUS	250.00	Uncon	
05/16/2005	Sentinal visit contracted in	878.11	Uncon	
12/15/2004	hypertension	401.1	Uncon	
12/15/2004	TYPE 2 DIABETES	250.00	Uncon	


Very little information usable for managing head injury!

(In spite of spending 100's of millions of dollars on the software.)


How do you make progress in managing the injured brain without good records (data)?

A New Medical Record for the Brain






A New Medical Record for the Brain




FDA Cleared
CE Mark




No FDA/CE

Current Medical Record




Brain Medical Record





Remote review of data



Data annotation tool



Switch Display: Monitor – Medical Record





Brain Monitor




Brain Medical Record

Apps





API can assign a region on any display for the app results

App Examples



- PRx
- Blood Pressure Management
- IBM Patient Priorities
- Prediction of Events
- Spreading Depolarizations

PRx – Loss of Autoregulation




Correlation
Coefficient

PRx


PRx – Loss of Autoregulation



PRx Correlation ABP-ICP assesses Cerebrovascular Reactivity (CVR)




Calculation of PRx




An Index of Cerebral Autoregulation

MAP





Average across 10 seconds of data to get one sample (non overlapping windows)




PRx

Compute Pearson's correlation coefficient across windows of 300 seconds, with an update rate of 60 seconds (80% overlap)

Mean ICP





One minute average at the cursor time.

Total number of minutes in the entire monitoring session.

Average value and number of minutes outside the desired goal, across the entire monitoring session.

The red area and % represent the relative amount of time spent outside the desired goal.

Metrics calculated within a shift.



Shifts are defined by going to **Tools** → **Settings...**, in the *Threshold Display* tab.

Alternatively, metrics can be calculated based on the display duration (by changing the **Time Base** in the measurement settings).

CNS Reader: Plug-in Architecture




Remote Monitoring of Blood Pressure




Courtesy – Dr. Greg Hawryluk

IBM – Prioritizing Patients




*** Connection Established ****


Spreading Depolarizations





Spreading depolarizations = class of pathologic waves characterized by near-complete sustained depolarization of neurons/astrocytes that propagate through gray matter at 1-5 mm/min





Physiological biomarker of brain injury.




SDs – Prolonged perfusion deficits



SDs - Correlates to other physiology



Preferred Platform




The preferred platform for slow potential EEG such as **Spreading Depolarization**

Recording



CNS
Advanced
ICU Amp




Electrocorticography

Display



CNS Monitor for Bedside Review
and Data Integration



Analysis




CNS Reader for Advanced Visualization

Easy Data Export for Advanced Analytics



SD App Display





Two Parallel Stories




1970s →



Brain Injury and Network Physiology




Brain Systems



What system interactions do we explore?


- Start simple
- Not all data on every patient
- Append to existing trial

Some ideas for Network Physiology



- NP can help quantify endpoint of TBI
 - Example is sleep after TBI
- Determination of Patient State
 - Detect a downward transition
- Origin of seizures
- Blood pressure management.
 - We lose the “network connection” ...the neurovascular coupling in TBI.

Conclusions



- Neuro is about 30 years behind cardiology
 - Due to complexity of the brain and monitoring technology
 - Look at what's going on in cardiac...this is where neuro will be
- We are just starting to get data that can answer questions
 - This area is ripe for investigation
 - We don't even know the questions to ask in some areas
- Funding seems to be available

Advice



- Starting a company, being an entrepreneur
 - Extremely rewarding, extremely risky
 - Never had a boss
 - Always worked for myself
 - Never been married
- My Mentor was/is Making Mistakes
- I am amazed at how much “luck” played in the success of my company
- One of the most important characteristics that helped me was “persistence”
- Follow your passion, the money will follow